



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁶ : G01N 31/22	A1	(11) International Publication Number: WO 95/13535 (43) International Publication Date: 18 May 1995 (18.05.95)
(21) International Application Number: PCT/US94/13416 (22) International Filing Date: 14 November 1994 (14.11.94) (30) Priority Data: 08/152,541 12 November 1993 (12.11.93) US (71) Applicant: BOEHRINGER MANNHEIM CORPORATION [US/US]; P.O. Box 50528, 9115 Hague Road, Indianapolis, IN 46250 (US). (72) Inventors: LEE, Debra, McDonald; 3654 E. Fallcreek Parkway N. Drive, Indianapolis, IN 46205 (US). BURKE, David, William; Apartment D, 1931 Madison Ct., Carmel, IN 46032 (US). TABB, David, Lee; 3144 W. 100 North, Indianapolis, IN 46140 (US). (74) Agents: AMICK, Marilyn, L. et al.; Boehringer Mannheim Corporation, 9115 Hague Road, Indianapolis, IN 46250 (US).		(81) Designated States: CA, JP, European patent (AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE). Published <i>With international search report.</i>
(54) Title: GLUCOSE CALIBRATOR AND CONTROL MATERIAL FOR TEST STRIPS (57) Abstract A non-serum based calibrator and control reagent is disclosed which is useful for calibrating and validating devices such as test strips for determining glucose. The reagent composition contains water, a predetermined amount of glucose, and a dihydroxy alcohol having more than 5 carbon atoms. An especially preferred dihydroxy alcohol is dipropylene glycol. A method of making the calibrator and control reagent is also disclosed.		

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AT	Austria	GB	United Kingdom	MR	Mauritania
AU	Australia	GE	Georgia	MW	Malawi
BB	Barbados	GN	Guinea	NE	Niger
BE	Belgium	GR	Greece	NL	Netherlands
BF	Burkina Faso	HU	Hungary	NO	Norway
BG	Bulgaria	IE	Ireland	NZ	New Zealand
BJ	Benin	IT	Italy	PL	Poland
BR	Brazil	JP	Japan	PT	Portugal
BY	Belarus	KE	Kenya	RO	Romania
CA	Canada	KG	Kyrgyzstan	RU	Russian Federation
CF	Central African Republic	KP	Democratic People's Republic of Korea	SD	Sudan
CG	Congo	KR	Republic of Korea	SE	Sweden
CH	Switzerland	KZ	Kazakhstan	SI	Slovenia
CI	Côte d'Ivoire	LI	Liechtenstein	SK	Slovakia
CM	Cameroon	LK	Sri Lanka	SN	Senegal
CN	China	LU	Luxembourg	TD	Chad
CS	Czechoslovakia	LV	Latvia	TG	Togo
CZ	Czech Republic	MC	Monaco	TJ	Tajikistan
DE	Germany	MD	Republic of Moldova	TT	Trinidad and Tobago
DK	Denmark	MG	Madagascar	UA	Ukraine
ES	Spain	ML	Mali	US	United States of America
FI	Finland	MN	Mongolia	UZ	Uzbekistan
FR	France			VN	Viet Nam
GA	Gabon				

GLUCOSE CALIBRATOR AND CONTROL MATERIAL FOR TEST STRIPSBACKGROUND

The present invention relates to calibrator and
5 control material useful in calibrating and validating
testing devices such as test strips and dipsticks. More
particularly, the present invention relates to a non-
serum based, aqueous glucose calibrator and control
material and to a method for making said calibrator and
10 control material.

The field of clinical chemistry and clinical
analysis is concerned, *inter alia*, with the determination
and quantification of various substances in body fluids.
Many examples of the substances which are determined can
15 be given, and these include cholesterol, urea, cations,
and glucose. These examples of analytes, as well as
others, are assayed in diverse body fluids such as urine
and blood.

The monitoring of the level of glucose in blood is
20 important to the management of diabetes. The level of
glucose in the blood is controlled by the amount of
carbohydrate ingested and by insulin. Too much insulin
lowers the glucose level, and too little will result in
an abnormally high level of glucose. Both circumstances
25 lead to serious health problems for the diabetic.

Most of the glucose testing done outside of the hospital laboratory is done in non-laboratory settings such as nurses' stations, physicians' offices and at home. Testing is frequently done by measuring the amount
5 of glucose in urine. As the level of glucose rises in the blood, it exceeds the ability of the kidney to reabsorb it, and glucose is excreted into the urine.

Although measurement of glucose in urine is useful, measurement of glucose in blood provides a more accurate
10 reflection of the condition of the subject. Urine glucose does not accurately reflect the level of glucose in the blood since the level of glucose in urine is determined by the level of glucose in the blood and the ability of the kidney to reabsorb the glucose.
15 Therefore, the urine sample cannot tell the diabetic how low his glucose level is.

Dry reagent test strips, sometimes referred to as dipsticks, are widely used for detecting glucose in urine and blood. These devices are characterized by their
20 simplicity of use. In general, such test strips comprise plastic strips provided at one end thereof with an absorbent paper portion which has been impregnated with reagents such as an enzyme system and a color indicator compound which produces or changes color to form a
25 detectable signal when the test strip is contacted with the analyte being determined. This change in color can

- be measured by comparing the color formed on the strip with a standard color chart calibrated to various glucose concentrations. More recently, however, to more accurately control the level of glucose in blood,
- 5 instruments have been developed which measure the color change in a reflectance photometer and thereby produce quantitative results. Examples of reaction systems which measure glucose using reflectance measurements include oxidative reactions, such as the glucose
- 10 oxidase/peroxidase method, and reductive reactions, such as the glucose oxidase/ferricyanide method. The latter method is described in detail in Freitag, U.S. Pat. No. 4,929,545, the content of which is herein incorporated by reference. Instruments have also been developed which
- 15 determine glucose by means of electrochemical methods in which a change in current is measured. An example of this technology is the amperometric, biosensor method described in PCT Application No. PCT/US90/07374, the content of which is herein incorporated by reference.
- 20 It will be understood that clinical analysis of the type described herein requires that any testing system be extremely accurate. In particular, when automated systems and instruments are used, it is essential to ensure that the elements of the analysis are reliable and
- 25 that the measurement taken is valid. It is for this purpose that calibrator and control reagents are used.

- 4 -

Westgard and Klee, in Textbook of Clinical Chemistry, N.W. Tietz, Ed., 1986, p. 430, define "reference material" as "a material or substance one or more properties of which are sufficiently well
5 established to be used for the calibration of an apparatus or for the verification of a measurement method." "Calibration and test material" is defined as "a reference material or solution with which the test sample is compared in order to determine the
10 concentration of analytes or other quantities."

The Westgard and Klee reference defines "control material" as "a specimen, or solution, which is analyzed solely for quality control purposes and is not used for calibration purposes." This standard reference work goes
15 on to describe some of the requisites of control materials as follows: "They need to be stable materials, available in aliquots or vials, that can be analyzed periodically over a long time. There should be little vial-to-vial variation so that differences between
20 repeated measurements can be attributed to the analytical method alone."

The above-cited reference, at page 433, discusses how the matrix of the control material should be the same as the material being analyzed. To that end, modified
25 human serum is discussed as one type of control material. Indeed, the art now recognizes the term "control serum"

- 5 -

as referring to control material based upon serum. This terminology will be used herein and is different from the term "control reagent," which, as used hereafter, refers to a control material which is not based upon, and which
5 does not contain, serum of any type.

As has been pointed out above, one of the criteria which control, and also calibration, materials have to satisfy is stability. Control materials based upon serum, however, are inherently unstable due to the
10 various components contained therein. Further, sera will vary from source to source, so uniformity from lot to lot cannot be guaranteed. Thus, it is sometimes desirable to have a calibrator or control material based upon a non-serum or serum-free medium.

15 An example of a serum-free control medium, or "control reagent" as used herein, is described in U.S. Pat. No. 4,729,959, issued to Ryan, which is directed to "a stable glucose reference control." This control contains glucose in a range of from about 40 to 500
20 mg/dl, together with fixed red blood cells, in an aqueous suspension. The range of glucose concentrations given are sufficient to cover just about all ranges of glucose found in, e.g., blood.

The Ryan '959 patent points to a problem with
25 aqueous control reagents at column 1, lines 50-55. Briefly, erythrocytes impart a degree of viscosity to

- 6 -

blood which is absent in water based systems. This problem was also recognized in U.S. Pat. No. 3,920,580 issued to Mast. This patent teaches that aqueous solutions had not been consistent, and that a lack of reproducibility was observed when dry reagent strips were used with such materials. Mast teaches that suitable reagents could be prepared by using an antidiffusing agent in combination with glucose and water. The antidiffusing agents taught by Mast include polyvinylpyrrolidone, polyvinyl alcohol, polyethylene glycol, dextran, and bovine serum albumin. Beneficial amounts are taught to be between about 3 and 35 percent of antidiffusing agent. The control solution may also include adjuvants to obtain a particular color or physical appearance, which include colored latex particles and water-insoluble lake dyes.

Maurukas, in U.S. Pat. No. 3,876,375, describes a reference control which is stable in the liquid state at temperatures ranging from -20°C to ambient room temperatures. Maurukas' compositions are prepared by removing from 20 to 40% of the water from an aqueous biological material and replacing it with substantially the same concentration of an alkylene polyol containing from 2 to 5 carbon atoms. Suggested alkylene polyols are ethylene glycol, propylene glycol, butylene glycol, pentanediol and glycerol. The compositions described are

- 7 -

said to be suitable for serum analysis using sequential multichannel automated analyzers. No mention is made as to the possible suitability of the compositions with test strips for measuring glucose using whole blood samples.

5 Kenamer *et al.*, in U.S. Pat. No. 5,028,542, the content of which is herein incorporated by reference, describe a non-serum based, glucose measurement control reagent in which the viscosity agent polystyrene sulfonate is used.

10 It has now been found that a suitable glucose calibrator and control reagent can be formed without using any of the materials referred to in the prior art as required ingredients. Rather, by combining a soluble dihydroxy alcohol with a predetermined amount of glucose
15 and water, along with additional optional materials, a suitable glucose calibrator or control reagent can be made.

SUMMARY OF THE INVENTION

20 The present invention is a non-serum based glucose calibrator or control reagent which comprises a predetermined, known amount of glucose, water, and a soluble dihydroxy alcohol containing more than 5 carbon atoms. A preferred dihydroxy alcohol is dipropylene
25 glycol, and a preferred concentration for the dihydroxy alcohol is between about 20 to 30 percent by weight of

the reagent composition. It was found, quite unexpectedly, that the composition of the present invention is useful in calibrating and validating testing devices such as test strips for the measurement of
5 glucose. Additional materials such as a buffer, a preservative, a biocide, an ionic salt, or a surfactant, either alone or in various additive combinations, may be mixed with the three required components. Another aspect of the present invention is a method for making the
10 calibrator or control reagent by mixing the glucose, water, and the dihydroxy alcohol containing more than 5 carbon atoms together.

Essential to the invention are a predetermined amount of glucose, water, and the recited dihydroxy
15 alcohol. The water is used, of course, to create a reagent solution in which the other components are dissolved. By "predetermined" is meant that, prior to formulation of the actual reagent, a concentration of glucose has been selected. This concentration may vary,
20 as those skilled in the art will recognize. As has been mentioned above, the art recognized, e.g., a range of from 40 to 500 mg/dl for control material, but one may envision broader ranges, e.g., from about 10 to 700 mg/dl, when the material is to be used for calibration
25 purposes.

- 9 -

The essential features of the invention, when the reagent is in the form of a solution, are the solvent (water), the predetermined amount of glucose, and the dihydroxy alcohol. The dihydroxy alcohol may be present
5 in, e.g., a range of about 20 to about 30 percent by weight of the reagent. The weight percent of the dihydroxy alcohol will be determined by criteria such as the final reagent viscosity desired and the desired diffusion or permeability characteristics of the
10 calibrator or control with the particular testing device with which it is to be used. Of course, the particular amount of dihydroxy alcohol selected should also be one which does not adversely interfere with the determination of glucose or have a negative effect on reagent
15 stability. It is not necessary that the control material have the same viscosity as whole blood; however, it is desirable that the permeability of the material, i.e., the diffusion rate of the glucose analyte, through the reagent matrix of the test strip approximate that of
20 whole blood.

Optional additional components of the calibrator or control reagent include typical additives such as buffers, preservatives, surfactants, and ionic salts. With respect to buffers, some preferred species are
25 succinate, HEPES (4-[2-hydroxyethyl-1-piperazine] ethane sulfonic acid), CHES (2-[N-cyclohexylamino] ethane

- sulfonic acid), MOPS (3-[N-morpholino] propane sulfonic acid), MEPS (2-[N-morpholino] ethane sulfonic acid), and CAPS (3-[cyclohexylamino-1-1-propane] sulfonic acid) buffers. Preferred preservatives or biocides include
- 5 imidazolidinyl urea, available under the trade name GERMALL 115 (GAF Chemicals Corp.), methylparaben (methyl-*p*-hydroxybenzoate) or methanol (((2-(dihydro-5-methyl-3(2H)-oxazolyl)-1-methylethoxy)methoxy)methoxy), available as COSAN 145 (Cosan Chemical Corp.),
- 10 phenoxyethanol and gentamycin sulfate, both individually and in combination.

- It may also be desirable to include a colored or colorable substance in the reagent mixture. This can be desirable because body fluid samples frequently possess a
- 15 particular color as one of their properties. As the control reagent is being used to calibrate per a body fluid sample, it can be useful to calibrate against conditions as similar to the tested fluid as possible, including color.

20

DESCRIPTION OF PREFERRED EMBODIMENTS

Example 1

- A preferred formulation of the calibrator and control reagent of the present invention was prepared
- 25 having the following composition:

HEPES, Na salt

30 mM

- 11 -

HEPES, free acid	30 mM
CaCl ₂	130 mM
NaCl	120 mM
GERMALL 115	0.3 %
Methylparaben	0.3 %
Water	77.5 %
Dipropylene glycol	22.5 %

This reagent was adjusted with NaOH to have a final pH of 7.5. Following preparation, the reagent was divided into aliquots and spiked with predetermined amounts of glucose. Theoretical or target glucose levels
5 selected were 10, 20, 40, 60, 80, 100, 130, 180, 300, 550 and 700 mg/dl. Each of the calibrator/control mixtures was then divided into vials, and duplicates were stored under a variety of temperature conditions.

Example 2

10 Each of the calibrator and control mixtures described in Example 1 was then tested for its efficacy after 1 week of storage at 25°C. As explained above, one of the most important features of a control reagent is its consistency, meaning that values obtained using it
15 should be fairly uniform from run to run.

With this in mind, each of the control materials stored at 25°C was applied to test strips containing the glucose determination system described in PCT/US90/07374. Briefly, this publication describes the determination of
20 glucose using an amperometric biosensor system.

- 12 -

Twelve replicates of each calibration/control mixture were measured, and the mean current readings at 10 seconds, standard deviations, and coefficients of variation were calculated. Reference glucose values for
5 the samples were determined using hexokinase methodology and an HITACHI 705 (Boehringer Mannheim Corp., Indianapolis, IN) analyzer. The results obtained are shown in Table 1 below.

Table 1

HK Glucose (mg/dl)	Mean (μ A)	Std. Dev.	Coeff. of Variation
11.0	1.30	0.06	4.75
41.5	3.81	0.07	1.84
63.0	5.85	0.09	1.60
85.0	7.90	0.10	1.31
106.0	9.98	0.12	1.17
137.5	13.18	0.21	1.63
191.0	18.30	0.16	0.87
316.5	30.29	0.92	3.04
573.5	52.40	2.31	4.41
711.5	64.05	3.81	5.96

These results show a level of consistency well within that required of a calibration and control reagent, as is indicated by the comparative standard deviation values and coefficients of variation reported for each set of tests.

It will be understood that the specification and examples are illustrative but not limitative of the present invention, and that other embodiments within the spirit and scope of the invention will suggest themselves to those skilled in the art.

What is claimed is:

1. A calibrator or control reagent for glucose determination comprising a mixture of a predetermined amount of glucose, water, and a dihydroxy alcohol having
5 more than 5 carbon atoms.
2. The calibrator or control reagent of claim 1, wherein said dihydroxy alcohol is dipropylene glycol.
3. The calibrator or control reagent of claim 1, wherein said dihydroxy alcohol is present in an amount
10 ranging from about 20 to about 30 percent by weight of said calibrator or control reagent.
4. The calibrator or control reagent of claim 1, further comprising a buffer.
5. The calibrator or control reagent of claim 1,
15 further comprising a preservative.
6. The calibrator or control reagent of claim 1, further comprising an ionic salt.
7. A calibrator or control reagent for glucose determination comprising a mixture of glucose, water,
20 dipropylene glycol, a buffer, a preservative, and an ionic salt.
8. A process for making a calibrator or control reagent for glucose determination comprising combining a predetermined amount of glucose, water, and a dihydroxy
25 alcohol having more than 5 carbon atoms.

9. The process of claim 8, wherein said dihydroxy alcohol is dipropylene glycol.

10. The process of claim 8, wherein said dihydroxy alcohol is present in an amount ranging from about 20 to
5 about 30 percent by weight of said calibrator or control reagent.

11. The process of claim 8, further comprising combining a material selected from the group consisting of a buffer, a preservative, and an inorganic salt with
10 said predetermined amount of glucose, water, and dihydroxy alcohol.

INTERNATIONAL SEARCH REPORT

 International application No.
 PCT/US94/13416

A. CLASSIFICATION OF SUBJECT MATTER

IPC(6) : G01N 31/22

US CL : 436/8, 12, 13, 14, 15, 16, 18, 95

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

U.S. : 436/8, 12, 13, 14, 15, 16, 18, 95

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practicable, search terms used)

STN Search terms: dipropylene glycol, control, calibrat7, glucose, dextrose

C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X --- Y	US, A, 3,920,580 (MAST) 18 November 1975, see entire document.	1,3,5-6,8,10 ----- 2,4,7,9,11
X --- Y	JP, A, 62-412 (SHIMIZU ET AL) 06 January 1987, see entire document.	1-2,5,8-9 ----- 7,11
Y	US, A, 5,028,542 (KENNAMER ET AL) 02 July 1991, see entire document.	2,4,7,9,11
Y	US, A, 3,876,375 (MAURUKAS) 08 April 1975, see entire document.	2,7,9,11
Y	JP, A, 3-112908 (URASHIMA ET AL) 14 May 1991, see entire document.	2,7,9,11

☒ Further documents are listed in the continuation of Box C.
 ☐ See patent family annex.

* Special categories of cited documents:	* later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
A document defining the general state of the art which is not considered to be of particular relevance	*X* document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
E earlier document published on or after the international filing date	*Y* document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art
L document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)	
O document referring to an oral disclosure, use, exhibition or other means	
P document published prior to the international filing date but later than the priority date claimed	*A* document member of the same patent family

Date of the actual completion of the international search

19 JANUARY 1995

Date of mailing of the international search report

08 FEB 1995

 Name and mailing address of the ISA/US
 Commissioner of Patents and Trademarks
 Box PCT
 Washington, D.C. 20231

Facsimile No. (703) 305-3230

Authorized officer

Arien Soderquist

Telephone No. (703) 308-0196

INTERNATIONAL SEARCH REPORT

International application No.
PCT/US94/13416

C (Continuation). DOCUMENTS CONSIDERED TO BE RELEVANT

Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
A	US, A, 4,729,959 (RYAN) 08 March 1988, see entire document.	1-11